






Endodontics / Endodontia / Endodoncia



ROCKMED® Industries **Endodontic Instruments** ROCKMED Endodontic Explorers have a long tip with a sharp point that is used to locate the opening of small canals during Endodontic procedures. ROCKMED Endodontic Explorers are made of spring tempered stainless steel so they retain their shape and resist breakage. MDN.E001 ROCKMED ( **MDN.E002** DG 16/23 ROCKMED CO 111111111111 11111111111 **MDN.E003** DG 16/17 ROCKMED CE REFERENCE talatalata talah **MDN.E004** 1 Abou-Rass Exp. A versatile explorer with two different, ROCKMED CE Listatatatata HILLIAM useful and sharp working ends. ROCKMED popular, basic apical explorer used to initiate the tracking groove in the root tip, and for general exploration and probing. **MDN.E005** Ø 0,7mm RI ROCKMED CE 11111111111 611111111 **MDN.E006 VXXXXXXXXXX** ROCKMED ( **MDN.E007** \*\*\*\*\*\*\* ROCKMED CERNANANA **ROCKMED Endodontic Excavators** have extra-long terminal shanks to reach deep into the cavity preparation. All are made of high quality stainless steel that is formed and precision ground by expert craftsmen, then hardened for the ultimate in sharp edge retention and durability. HILLIAM ROCKMED CE 11111111111 **MDN.E008** Excavator: The shank is long to reach canals. To curettage inside of tooth to base of pulp chamber. Sparkle 2.0 1.0 MDN.E011 **MDN.E009** 20 31W Ovel Spoon **MDN.E012 MDN.E010 MDN.E013** ROCKMED CE 2.0 **MDN.E014** ROCKMED CEN WWW.ROCKMED.CO Page 2









# **Endodontic Instruments**






Pluggers: are used to condense filling materials into cavity preparation. Ends are serrated. Burnishers: Used to condense, smooth and polish amalgam.







# ROCKMED Industries



# **Endodontic Instruments**

### MDN.E135

DOVGAN MTA CARRIER #DOVGAN16, Ø1.6MM, BENDABLE TIP, COMPONENT BASED



# RI ROCKMED

Simple Socket MDN.E136S

RI ROCKMED



R ROCKMED CE



MIRROR HANDLE, WITH ENDODONTIC FILE MEASURING SLIDE

Simple Socket MDN.E137S Cone Socket MDN.E137C

Cone Socket MDN.E136C



B) POCKMED CE

Serrated Tip For Reduced Slippage



# RI ROCKMED

MDN.200.01R
Carr Endodontic Retractor 90° Angled

ROCKMED

RI ROCKMED CE

### MDN.200.02R

Carr Endodontic Retractor 45° Angled Serrated Tip For Reduced Slippage



## **Endodontic Aspirators**

Modified Frazier style aspirator.

Tapered tip. Short dental style bend.

RI ROCKMED

MDA.ASP.24 1.5mm MDA.ASP.25 2.0mm MDA.ASP.26 2.5mm MDA.ASP.27 3.0mm MDA.ASP.28 4.0mm



### MDA.ASP.11 Endodontic aspirator made of 1/8" diameter stainless steel. 1.7mm opening

ROCKMED



Micro suction for root canal and root tip surgery. Fits standard high volume cut off valve (11mm). Endodontic aspirator made of 1/8" diameter stainless steel.


MDA.ASP.12 1.7MM ENDO ASPIRATOR MDA.ASP.13 2.3MM ENDO ASPIRATOR

# **Endodontic Forceps**





Micro forceps for removing broken files 75 mm



- The lid has a locking system. Easy lock-on cover secures instruments in container.
- It can store 48 files / instruments (classical instrumentation or nickel-titanium).
- It resists all modes of decontamination.
- It can be placed directly in the ultrasonic box.
- The preservation of the sterilization is allowed by the bagging of the endo box.
- Holds 48 Files




Teflon, Holds 48 Files







# **Root Canal Instruments**




Page12 WWW.ROCKMED.CO






# **Root Canal Instruments**



The Bone Graft Kit includes a number of essential dental instruments to be used in the placement of allograft to the surgical site.

There is no more accurate way to ensure replicable packing and placement of your bone grafting materials. This kit includes packers and carriers with serrated titanium tips:



# Guarantee: All ROCKMED instruments are guaranteed to be free of manufacturing and workmanship defect. Any instruments found to be defective, excluding defects due to misuse or mishandling, shell be repired or replaced free of cost. We recommend the use of proper and standard instrument cleaning, handling and preserving techniques to increase the life and functionality of delicate dental instruments.


The ROCKMED Logo, and the Name are Registered and Trade Mark of ROCKMED Industries



Scan our mobile website to your smartphone for quick reference



Corporate Head 61-B Liaqat Ali Road, Industrial Estate, Sialkot, Pakistan Phone # +92-523-555526





+ 92-333-8623458

492 331 3177777

info@rockmed.co

www.rockmed.co

# Precision Craftsmanship

**CERTIFIED COMPANY** ISO-9001 ISO-13485



Copyright © ROCKMED Industires Reprinting of this catalog even partially is not allowed. ROCKMED Management Information System.