


Skin Graft Mesher Skin expansion system

ROCKMED Skin Graft Mesher is a surgical skin graft tool designed for simple operation. Durable, stainless-steel construction for years of consistent performance. Simple to change-out cutter options. Easy ratchet handle motion for smooth transport of graft carrier through the mesher.



Mesher for skin grafting

Skin Graft Mesher is a tool for preparing skin for grafts.

Patches of skin that were previously collected with a dermatome are stretched with a mesher over a dedicated mesh, allowing the skin to be stretched in width. The reinforced surface of the transplanted skin spares the patient from further procedures. Vacant spaces in the transplanted part of the skin overgrow during the healing process.



