


Magnetic Drape / Mat



ROCKMED
Industries



Magnetic Drape / Mat

Code# **SIL2.MAG** Magnetic Mat 12" x 15" / 30 cm x 38 cm

Our Magnetic Mat provides a sterile, neutral field for the safe placement and hands-free transfer of instruments. It is latex free, reusable, and autoclavable (up to 134° C).

*Conforms to any surface shape given its flexible and non-slip backing drape
Round magnets to avoid rubber perforations*

Raised magnets hold instruments firmly in place while allowing for easy grasping


PRECISION CRAFTSMANSHIP

ISO 9001 2015
ISO 13485 2016



www.rockmed.co





ROCKMED Silicone Pin Mat

Silicone Pin Mat - SILICONE TRAYS INSERT MATS

- ~ Reusable - Autoclavable - Latex free
- ~ Highly flexible
- ~ Easy to cut
- ~ Embossed on the underside for optimum air circulation in the autoclave.
- ~ **FRONT SIDE** Pins keep instruments securely in place
- ~ **BACK SIDE** Numerous holes and elevated nubs provide superior venting and reduce water collection

GOOD TO KNOW: CUSTOM SIZES ARE AVAILABLE UPON REQUEST

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization.

The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.



Codes	Description	Sizes mm
SIL2-2025	Embossed silicone mat	197 x 197mm
SIL2-1405	Embossed silicone mat	198 x 86mm
SIL2-3215	Embossed silicone mat	337 x 210mm
SIL2-1055	Embossed silicone mat	445 x 191mm
SIL2-1215	Embossed silicone mat	465 x 212mm
SIL2-1914	Embossed silicone mat (without holes)	483 x 356mm
SIL2-4005	Embossed silicone mat	502 x 236mm
SIL2-6005	Embossed silicone mat (2 units)	536 x 249mm

